

Adapter reference

The step-by-step guide is [VENDOR-INTEGRATION-GUIDE.md](#). This page is the short reference.

What an adapter may import

These are kept stable:

Module	What you use
<code>openbedside.adapters.stream</code>	<code>StreamAdapter</code> : subclass it, write <code>connect</code> , <code>read_message</code> , <code>normalize</code> , <code>close</code>
<code>openbedside.adapters.base</code>	<code>DeviceAdapter</code> , if you need full control of <code>run()</code>
<code>openbedside.model</code>	<code>Device</code> , <code>Module</code> , <code>Channel</code> , <code>InfusionSource</code> , <code>Quantity</code> , <code>PumpStatus</code> , <code>SourceRole</code> , <code>utcnow</code>
<code>openbedside.validate</code>	<code>check_device(dev)</code> for your tests; levels <code>ERROR</code> , <code>WARNING</code> , <code>INFO</code>

The snapshot

```
Device      device_id, vendor, model, observed_at, device_clock, patient_weight,
            reported_patient_id (only an id the device itself holds), kind
├─ Module   module_id                (one per pump module or channel)
├─ Channel  channel_id, status, actual_rate (exactly one per module)
├─ InfusionSource role PRIMARY or SECONDARY, drug_name, concentration,
                  rate, dose_rate, vtbi, vtbi_remaining, volume_delivered
```

`Quantity(value, ucum)` . Supported UCUM units are the keys of `codes.UNITS` : `mL` , `mL/h` , `mg` , `mg/mL` , `ug/kg/min` , `[iU]/mL` , `[iU]/h` , `kg` . Rates are mL/h and volumes are mL.

`PumpStatus` : `INFUSING` , `KVO` , `PAUSED` , `IDLE` , `COMPLETE` , `ALARM` , `UNKNOWN` .

What an adapter must not set

`Device.patient` , `ehr_id` , `eui64` and `location` are set by the gateway from the device registry and the association step. `openbedside check` reports an `ERROR` if an adapter sets `patient` . A device the hospital has not registered is held, not sent.

Registering an adapter

```
[gateway]
adapter_paths = ["."]          # folders to search, relative to this config file

[[adapters]]
type = "mymodule:MyAdapter"    # module:ClassName
# any other keys are yours, available as self.config["..."]
```

Programming

Real-device adapters do not receive programming requests in this release. Leave `supports_programming` `False` and set `[orders] enabled = false`.